

Store Section

Dr. BRA IRCH

Ref No: F. No. 23/IRCH/SSK/L.O/16-17

Date: 17/05/2016

Subject: Purchase of Channel Pipettes, Mfg.by eppendorf AG Germany for Lab.Onco. Dept. (Dr. BRAIRCH)

The request letter received from Dr. Anita Chopra duly forwarded by Dr. Rajive Kumar (Prof. & Head Lab Oncology) for the procurement of Channel Pipettes, on proprietary basis from M/s eppendorf AG Germany

Sr. No.	Description Item /Cat No#
1.	3122 000.019 Research Plus, 8-channel 0,5-10µl Fully Autoclavable Spring Loaded Tip Cone: Improved ergonomics, Tight fit to the tip (not in 5ml and 10ml pipettes) Adjustment opening: Adjust your pipettes to a specific liquid and volume. Control Button: Very low operating force, Color indicates pipette volume, Positioned for perfect ergonomics. Removable independent channels: Gain flexibility for individual needs Ejector: Very low operating force, positioned for perfect ergonomics. Volume Display: 4 Digits, Magnifying shape. Eppendorf Perfect Piston System: Ultra light system made of Forton. Quick Connection Clip: Remove lower part easily
2.	3122 000 0.35 Research Plus , 8-channel 10-100µl Fully Autoclavable Spring Loaded Tip Cone: Improved ergonomics, Tight fit to the tip (not in 5ml and 10ml pipettes) Adjustment opening: Adjust your pipettes to a specific liquid and volume. Control Button: Very low operating force, Color indicates pipette volume, Positioned for perfect ergonomics. Removable independent channels: Gain flexibility for individual needs Ejector: Very low operating force, positioned for perfect ergonomics. Volume Display: 4 Digits, Magnifying shape. Eppendorf Perfect Piston System: Ultra light system made of Forton. Quick Connection Clip: Remove lower part easily.
3.	3881 000.031 PCR-Cooler 0,2 ml Blue Application For the protection, transport and storage of your sensitive samples and of course for the safe stoppage of reactions. Product features Cools sensitive samples reliably for more than an hour at 0°C (with 2 hour precooling at -20°) Temperature safety through color change when temperature is exceeded: > 7 °C : pink or light blue < 7 °C : violet or dark blue Simple, space saving Instead of ice bath: safe from contamination No freezing of samples with special cone contour Each individual bore shows the temperature of each single tube In the 96 - well format For 0.2 ml and 0.5 ml, tube strips, plates and single tubes LESS

The user department has submitted proprietary Article certificate in addition to one submitted by company.

The above documents are being uploaded on AIIMS WEBSITE for open information to invite objection, comments, if any, from any manufacturer regarding proprietary nature of the item within 15 days from the date of issue/uploading of the notification giving reference F. No 23/IRCH/SSK/L.O/16-17. The comments should be sent to stores officer, Dr. BRAIRCH AIIMS on or before 03/06/2016, up to 4:30 p.m failing which it will be presumed that any other vendor/company is having no comment to offer and case will be decided purely on merits.


Stores Officer
Dr. BRA IRCH

Encl:

- PAC Certificate (enclosed)

Date: 22nd October, 2012

TO WHOM IT MAY CONCERN

We here by certify that "Eppendorf Research Plus Pipette" is the proprietary in nature because of its specification & Special features like (Spring loaded TIP Core; Adjustment window to handle viscous solutions; Fortron Piston) manufactured by M/s Eppendorf AG Germany.

To the best of our knowledge; no other manufacturer world wide has similar/ identical Specification.

Yours faithfully,

Dr. Andrés Jarrin
Area Sales Manager
Eppendorf AG
Barkhausenweg 1
22339 Hamburg
Germany

eppendorf^{AG}
Global Sales
Barkhausenweg 1
22339 Hamburg
Germany

Phone: +49 40-5 38 01-172
Fax: +49 40 5 38 01-840
jarrin.a@eppendorf.de



◀ Vorstand:
Dr. Dirk H. Ehlers (Vorsitzender)
Dipl.-Kfm., Detmar Annemermann
Dr. Heinz Gerhard Köhn
Dr. Michael Schroeder

Vorsitzender des Aufsichtsrats:
Dipl.-WirtschaftsIng. Klaus Fink

◀ Eppendorf AG
Barkhausenweg 1, 22339 Hamburg
Postanschrift: 22331 Hamburg
Tel.: (040) 5 38 01-0
Fax: (040) 5 38 01-556
www.eppendorf.com

Amtsgericht: Hamburg HRB 76249

◀ Deutsche Bank AG
Kto. 0202200 / BLZ 200 700 00
IBAN DE51 2007 0000 0020 3000 0000
S.W.I.F.T. DEUT DE 33

HSBC Trinkaus & Burkhardt AG
Kto. 800547003 / BLZ 300 000 000
IBAN DE56 3003 0500 0000 0000 0000
S.W.I.F.T. TUBD DE 33

Dr. Anita Chopra, MD, DM
Professor
Laboratory Unit
Dr. BRAIRCH, AIIMS, New Delhi